

Gynaecology & Proctology

NO#	Components	Shoulder capacity(ml)	Brimful capacity(ml)	Height (mm)	L*W (mm)	Weight (g)	Materials	Neck size (mm)	Total Height(mm)	Item code
1	Applicator	50	55	102.7	23.0*23.0	4.5	PP	20	144	BONA5775A6E
	Bottle			71.8	43.5*43.5	6.5	LDPE			
2	Applicator	100	114	102.7	23.0*23.0	4.5	PP	20	185.5	BONA5539ABE
	Bottle			97.5	52.0*52.0	15.4	LDPE			
3	Applicator	500	520	108.2	34.4*34.4	9.5	PP	30	233.2	BONA5530A6V
	Bottle			146.0	88.0*88.0	62.0	LDPE			
4	Applicator	150	160	79.2	22.5*22.5	2.5	PE	20	190	BONA4R0271F
	Bottle			121.78	50.0*50.0	16.0	LDPE			

NO#	Dose form	Capacity(g)	Length(mm)	Diameter(mm)	Materials	Item code
5	Gel	3	84.37	16.0*27.8	PP	BONA4G0023F
6	Gel	5	94.0	16.0*27.8	PP	BONA4G0006F
7	Gel	10	102.0	21.65*31.55	PP	BONA4G0026F
8	Gel	5	97.2	19.1	LDPE	BONA4G0045F

Main Parameters:

Neck size: 20mm, 30mm

Douche bottle volume: 50ml, 100ml, 150ml, 500ml

Gel applicator & tube: 3g, 5g, 10g

Materials: HDPE, LDPE, PP

Key features:

Comfortable design with smooth applicator surface to avoid self-injury

Compatible with gel and cream

Non-sensitizing, non-toxic and non-irritating

Vaginal applicator can be prefilled

Comply with CP, EP, USP food & drug direct contact regulations

